



**DV-60**

LED Light Curing Device

---

OPERATION    MANUAL

URIT Medical Electronic Co.,Ltd.



# CONTENT

Copyright and Notice .....	2
1. Introduction .....	3
2. Structure .....	4
3. Technical Specifications .....	6
4. Installation and disassembly .....	8
5. Operating Instructions .....	10
6. Precautions .....	13
7. Contraindications and applicable groups .....	15
8. Maintenance .....	15
9. Troubleshooting .....	16
10. Storage and Transportation .....	18
11. After-sale Service .....	18
12. Environmental Protection .....	18
13. Symbol Instructionsn .....	19
14. Statement .....	20
15. Attachment Lists .....	20
16. EMC .....	20

## **Copyright and Notice**

URIT Medical Electronic Co.,Ltd. ©All rights reserved.

Congratulations on becoming a respected customer of URIT Medical Electronic Co.,Ltd. and welcome to use DV-60 LED Curing Light, it will bring you a new experience and convenience.

This manual is compiled according to the relevant laws and regulations of China and the specific conditions of the DV-60 LED Curing Light manufactured by URIT Medical Electronic Co.,Ltd.

This manual contains the latest information as of the printing of this manual. URIT Medical Electronic Co.,Ltd. is solely responsible for the revision and explanation of the simplified English version of this manual, and reserves the right to change the relevant content after printing this manual without prior notice.

The pictures involved in this manual are schematic diagrams and are for reference only. If the pictures do not match the actual objects, the actual objects shall prevail.

All materials in this manual are protected by copyright law. Without the prior written consent of URIT Medical Electronic Co.,Ltd.,any content in the manual shall not be reproduced, copied or translated into other languages in any form.

The operator must operate in strict accordance with the operating instructions of this manual. Otherwise, URIT Medical Electronic Co.,Ltd. is not responsible for any errors and instrument



failures caused by illegal operations.

**Notice: URIT Medical Electronic Co.,Ltd. does not promise that the instrument is used for a particular purpose and makes any implied warranties for its merchantability and applicability.**

If you need after-sales service support, please contact Guilin Veirun Medical Technology Co., Ltd. or an authorized agent.

## 1. Introduction

### 1.1. Principle of operation

The device illuminates photosensitive resin to solidify it quickly based on the principle of light radiation.

### 1.2. Intended use

The device is intended to be used on a variety of dental compositions within the oral cavity.

### 1.3. The Features of Products.


- Stable constant light output power so that the solidification effect is not affected by the

consumption of remaining power.

- Parallel light ensures more focused output energy and better solidification effect.
- Efficient photovoltaic conversion. Once fully charged, it can be used for more than 300 times continuously under 10s working of Bright Light Mode, from full charge to low-battery alert.
- Ergonomics design, excellent and comfortable handle feel during the treatment.
- Aluminum alloy body, resistance to shatter and drop.
- Light source featuring 360° rotation to realize curing from all directions.
- Replaceable light source and enables detection of dental caries.
- The charge cradle is made of aluminum alloy to increase stability.
- Solidified lens fixed on light source, which are bonded with all-ceramic veneer and adhered to crown and fiber-reinforced composite root canal post to make solidification effect better.
- Replaceable Lithium battery, Auto power-off

## 2. Structure

LED Light Curing Device consists of main unit, charger cradle, battery, Lamp guide, shading plate and power adapter.



### 3. Technical Specifications

3.1. Dimensions:  $\phi$  18 mm x 228 mm.

3.2. Main unit net weight: 90 g.

3.3. Configurations: see the packing list for details,.

3.4. Power supply:

- Power supply classification:

Power supply by rechargeable lithium battery.

- Rechargeable lithium battery:

Standard voltage:DC3.7 V, capacity:800 mAh, with protection for overvoltage, overcurrent and short-circuits.

- Power Adapter:

Input: 100 V–240 V~, 50/60 Hz, 0.5A Max

Output: DC 5V  $\equiv$  1A

Fuse: T1AL 250 V

### 3.5. Characteristics of LED lamp:

10 W high–power blue LED lamp.

Lamp classification group: Risk Group 2.

Wavelength: 385 nm~515 nm.

Inspection method: When the LED lamp emits during normal operation,the lamp is in good condition.

Radiation during continuous working:  $\geq 400 \text{ mW/cm}^2$ .

Optical effective area of light source:  $75 \text{ mm}^2$ .

LED Light Curing Device is applicable to some common dental resinous materials in dental clinic,such as 3 M, Dentsply etc.

### 3.6. Working environment:

Environment temperature:  $5 \text{ }^{\circ}\text{C} \sim 40^{\circ}\text{C}$ .

Relative humidity:  $\leq 80\%$ .

Atmosphere temperature: 75 kPa~106 kPa.

### 3.7 Equipment Safety classification:

Operation mode: Short-time operation equipment.

Type of protection against electric shock: Class II equipment.

Degree of protection against electric shock: Type B application part.

Degree of protection against ingress of water: Ordinary equipment(IPX0).

Degree of safety of application in the protection of a Flammable Anesthetic Mixture with air,Oxygen or Nitrous Oxide: Non-AP,APG type equipment.

#### 4. Installation and disassembly

4.1. The top of the body unit featuring 360° rotation is replaceable.

As shown in figure.




4.2. Press and install the shading plate on the light guiding head.

4.3. Battery Replacement

Holding the front of the host and the battery compartment, screw it off counterclockwise, then remove the old battery from the battery, change a new battery, and then replace the battery storehouse on the host. The battery positive and negative can be normal used, no security risks.As

shown in figure:



4.4 When device needs to charge, firstly insert the DC connector of the power adapter into the charging jack marked on the charger cradle, and then connect the power adapter to the AC 100V–240 V electricity network, and then place the host in the center of the charger cradle to charge.

## 5. Operating Instructions

### 5.1. Orthodontic mode

Press button "M/T" for 1 second and release the button when the buzzer warns for once. Then you can choose 3s or 5s working time. The working cycle of orthodontics is 10 times, and the screen displays "x10", and the LED lamp is automatically turned off after 10 times of cyclic irradiation.

5.2. Press button "M/T" for 2 second and release the button when the buzzer warns for twice.

The four output of power density can be chosen as following:


- High Power Density: The screen shows "P1", The power density is about  $1800\text{mW/cm}^2 \sim 2100\text{mW/cm}^2$ .
- Ultra High Power Density: The screen shows "P2", The power density is about  $2000\text{mW/cm}^2 \sim 2300\text{mW/cm}^2$ .
- Standard Power Density: The screen shows "P3", The power density is about  $1000\text{mW/cm}^2 \sim 1300\text{mW/cm}^2$ .
- Caries Detection Power Density: The screen shows "P4" and "CHECK", and the purple light is output independently. The purple light works under proper use can be regarded as normal.



5.3. Tap button "M/T" to select working time interval:

- In "P3" mode, the time interval can be chosen from 1,2,3,4,5,10,15,20,25,30, 35, 40 second.
- In "P2" mode, the time interval can be chosen from 1,2,3,4,5,10,15,20 second.
- In "P1" mode, tap button "M/T" without any action. The time is set to 3s or 5s, and the operation method for switching is the same as 5.1.
- In "P4" mode, the time interval can not be chosen.

Information on the screen is shown in the figure:



5.4. During the operation, put the disposable sleeve on the top of the main unit, aim the top at the correct position, press the ON/OFF button ( "P" )and the main unit will produce "Di"sound, the curing light radiates blue light and starts working according to the set modes. Meanwhile, it starts counting down from the set working time interval, it stops working when counting down to "0".The screen displays the setting working time interval again.

5.5. Operation can be stopped by press the ON/OFF button ("P").

5.6. After a working cycle, operator can press the ON/OFF button ("P") to start another working cycle. Stop operating if the equipment is burning obviously, let equipment cool down before restarting. Suggest continuous working cycle less than 5 times.

5.7. Battery Indicator: Low power detective circuit is fixed inside of the main unit.If battery indicator in the screen has only 1 bar left, please charge in time.

5.8. When the battery needs to be charged, insert the DC plug of the power cable into the charging jack marked by the charge cradle, and then connect the power adapter to the AC 100V-240V power supply. At this time, the indicator light of the charge cradle is green, indicating that it enters the standby state. After placing the host in the center of the charging station, the charging station indicator lights up orange, indicating that charging

begins; When charging finished, the indicator light at the bottom of charger turn green.

5.8. After operating, take off the disposable sleeve and throw away, avoid reusing. Power output will be decreased by 5–10% if using the disposable sleeve.

5.9. The product will turn off automatically if there is no operations within 2 minutes. Turn it on by pressing any button.

5.10. The depth of solidification of dental resin composite is no less than 4mm per 10 second. The recommended separation distance between luminous point and solidifying point is 2mm.

## 6. Precautions

**Warning: It's our duty to provide users correct usage rules and safety notices.**

6.1 Using the product according to the instruction manual, other tasks are out of support.

6.2 Please charge the battery at least 4 hours before first time usage.

6.3 In order to prevent cross-infection, it is forbidden to reuse the disposable sleeve.

6.4 The top of the main unit can be turned by 360° and it is demountable and replaceable.

6.5 DV-60 LED light curing device only used by professionally trained people.

6.6 Please put it away from the touch of kids.

6.7 During operation, the light should be aimed straightly at the dental resin composite to ensure the effect of solidification.

6.8 Avoid aiming the light at eyes, please use brake sack or protective goggles supplied by our company to protect your eyes.

6.9 Please use the power adapter which is designed and supplied by our company. It may cause potential dangers to lithium and control circuit by using the power adapter designed or supplied by other manufacturers.

6.10 It is forbidden to put metal on the charger because it may burn the internal circuit. Unplug the plug of charger when not charging.

6.11 Please charge the battery in cool and ventilated room.

6.12 The product should not be used more than rated working time set in the instruction manual to avoid damaging teeth for the high temperature. Shut down and cool the product when it is used for 5 times continuously.

6.13 It is forbidden to extrude, shake or rock the battery. It is forbidden to self-taking apart the battery, in order not to result in short-circuit or leakage and it is forbidden to put the battery with metal.

6.14 Long time do not use the product, please take the lithium battery apart from the main unit. It is recommended that the product should be installed at the medical site where there are no high frequency high voltage equipment within 5 meters to ensure the product work properly.

## **7. Contraindications and applicable groups**

7.1 People suffered ophthalmic surgeries or sensitive to light, pregnant women, children and the heart disease patients should not use the product.

7.2 Patients with retinal diseases should be cautious to use the product.

7.3 Patients with photosensitization and solar dermatitis or using photosensitive drugs should not use the product.

Applicable to all groups except contraindications.

## **8. Maintenance**

8.1 This product does not consist of the self-maintainable spare parts. The maintenance of this product should be taken by the appointed professional or special repair shop.

8.2 Users can change brake sack, light source and lithium battery. Please use accessory which is designed and supplied by our company. Contact with the local dealer or our company if you want to buy. It may cause potential dangers to curing light or other damages which is designed and supplied by other manufacturers.

8.3 The accessory of the product should be cleaned by clean water or neutral sterilized liquid. Do not soak. Do not use highly volatile and different solvent to clean this product, which can cause the signs on the control pane to fade.

8.4 Please check whether there are any remains on luminous point and clean the top of the main unit with 75% alcohol tampon after using to avoid pollute and ensure solidified effect.

8.5 Please charge in time when battery level is low; get lithium battery out for safe-keeping when not using the product for a long time.

## 9. Troubleshooting

Faulty	Possible cause	Solutions
No indication, no response	1.Battery is out of power. 2.Faulty of battery. 3.The main unit battery	1.Charge the product/change battery 2.Change battery

	protection system works	3.Place the main unit on the charge for activation
"Er"shown on the screen.	Faulty of main unit.	Send to after-sale service for repair.
Light intensity is weak.	There is resin on the top of the main unit.	Clear the resin.
The equipment is not charging when the adapter is connected.	1.The adapter is not connected well. 2.Faulty from adapter or incompatible. 3. The charger is out of order.	1.Reconnect. 2.Change the adapter 3.Send to after-sale service for repair.
Usage time shortened on a single full charge.	Smaller battery capacity.	Change battery.

**Note: If such solutions are completed, the product still cannot work normally, please contact with the distributor or our company.**

## 10. Storage and Transportation

10.1 The product should be handled carefully and lightly. Be sure that it is far from the vibration, and installed or kept in a cool, dry and ventilated place.

10.2 Do not store the product together with the combustible, poisonous, caustic or explosive goods.

10.3 The product shall be stored at the location as follow:

- ① Relative humidity:  $\leq 90\%$ , ② Atmospheric pressure: 75 kPa~106 kPa
- ③ Temperature:  $-20\text{ }^{\circ}\text{C}$ ~ $55\text{ }^{\circ}\text{C}$ .

10.4 During transportation, excessive shock and vibration shall be prevented, and handled carefully and lightly. Keep it away from sunlight, snow or rain. It shall not be packed with dangerous goods.

## 11. After-sale Service
































From the date this product has been sold, base on the warranty card, we will repair this equipment free of charge if it has quality problems, please refer to the warranty card for the warranty period.

## 12. Environmental Protection

There is not any harmful element in our product. It can be disposed according to the local law.



## 13. Symbol Instructions

 Manufacturer's logo	 Type B applied part	 Refer to Instruction manual/ booklet	 Caries Detection Mode
 Environment-friendly use period	 Manufacturer	 Use-by date	 Standard Power Mode
 Serial number	 Do not dispose of the product into the ordinary municipal waste or garbage system	 Caution	 Ultra High Power Mode
 Atmospheric pressure limitation	 Fragile, handle with care	 Rotating plug	 High Power Mode
 Humidity limitation	 Temperature limit	 This way up	 Medical Device
 Keep away from rain	 Class II equipment	 Direct current	 Alternating current
 Full Power Mode	 Progressive Mode	 Pulsed Mode	 Battery Indicator
 Date of manufacture	 Mode/Time Button	 On/Off Button	

## 14. Statement

All rights of modifying industrial design, inner structure, instruction manual, etc, of the product are reserved to the manufacturer without further notice. The pictures are only for reference. The final interpretation rights belong to URIT Medical Electronic Co.,Ltd..

## 15. Attachment Lists

No.	Attachment Name	Pcs
1	Main unit	1
2	Lamp guide	1
3	Power adapter with one USB cable	1
4	Charger cradle	1
5	Rechargeable lithium battery	1
6	Light source	1

## 16. EMC

**Note:**

**(1) Unauthorized changes or modifications without the affirmative consent of URIT Medical Electronic Co.,Ltd. may cause EMC problems to the product or other equipment.**

**(2) The DV-60 LED light curing device has been tested in accordance with operating procedures of EMC.**

#### 16.1 Requirements of cable installation

Cable Name	Cable Type	Cable Length
Power supply output line.	Unshielded parallel line.	1.5 meter.

#### 16.2 Key parts of EMC

Key parts of EMC of the product are LED driver chip and the power adapter. Using or replacing accessories which are not designed and supplied by our company would result in performance degradation of electromagnetic emissions and electromagnetic immunity. Therefore, do not replace parts of the product without permission.

#### 16.3 Guidance and Manufacturer's Declaration--Electromagnetic Emissions.

Guidance and Manufacturer's Declaration--Electromagnetic Emissions.

The DV-60 LED light curing device is indicated for use in the electromagnetic environment specified below. The customer or the use of the models DV-60 LED light curing device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment --guidance.
RF emissions GB 4824	Group 1	The models DV-60 LED light curing device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions GB 4824	Group B	The models DV-60 LED light curing device are suitable for used in domestic establishment and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

Harmonic emissions GB 17625.1	Not applicable	
Voltage fluctuations /flicker emissions GB 17625.2	Applicable	

#### 16.4 Guide and Manufacturer's Statement--Electromagnetic Immunity.

Guide and Manufacturer's Statement--Electromagnetic Immunity			
The DV-60 LED light curing device is indicated for use in the electro-magnetic environment specified below. The customer or the use of the models DV-60 LED light curing device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment--Guidance.
Electrostatic discharge(ESD) GB/T 17626.2	$\pm 6$ kV contact discharge $\pm 8$ kV air discharge	$\pm 6$ kV contact discharge $\pm 8$ kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Electrical fast transient/burst GB/T 17626.4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge GB/T17626.5	$\pm 1$ kV line to line $\pm 2$ kV line to earth	$\pm 1$ kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines GB/T 17626.11	<5% UT ( > 95% dip in UT) for 0.5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT (95% dip in UT) for 5 sec.	<5% UT ( > 95% dip in UT) for 0.5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT (95% dip in UT) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the models DV-60 LED light curing device require continued operation during power mains interruptions, it is recom

			mended that the models DV-60 LED light curing device be powered from an nonnot rruptible powersupply or a battery.
Power frequency (50/60 Hz) magnetic field GB/T17626.8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note:UT is the a.c. Mains voltage prior to application of the test level.			

## 16.5 Guide and Manufacturer's Statement--Electromagnetic Immunity.

## Guide and Manufacturer's Statement--Electromagnetic Immunity

The DV-60 LED light curing device is indicated for use in the electro-magnetic environment specified below. The customer or the use of the models DV-60 LED light curing device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment--Guidance
Conducted RF GB/T17626.6	3 Vrms 150 kHz~ 80 MHz	3 Vrms	Portable and mobile RF Communications equipment should be used no closer to any part of the models DV-60 LED light curing device, including cables, than the recommended separation distance calculated from the frequency of the transmitter. Recommended Separation Distance.
Radiated RF GB/T17626.3	3 V/m 80 MHz~2.5 GHz	3 V/m	



			$d = \left[ \frac{3.5}{V_1} \right] \sqrt{p}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{p} \quad 80\text{MHz} \sim 800\text{MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{p} \quad 800\text{MHz} \sim 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter In watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m).</p> <p>Field strengths from fixed RF transmitters, as determined by</p>
--	--	--	--

			<p>an electromagnetic site survey,</p> <p>a. should be less than the compliance level in each frequency range.</p> <p>b. Interference may occur in the vicinity of equipment marked with the following symbol:</p>
--	--	--	--



Note1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note2: These guidelines may not apply in all situations. Electromagnetic propagation is

affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the models DV-60 LED light curing device are used exceeds the applicable RF compliance over, the model DV-60 LED light curing device level should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the models DV-60 LED light curing device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

16.6 Recommended separation distances between portable and mobile RF communications equipment and the models DV-60 LED light curing device.

Recommended separation distances between portable and mobile RF communications equipment and the models DV-60 LED light curing device.

The models DV-60 LED are intended for use in electromagnetic environment in which radiated RF disturbances is controlled. The customer or the user of the models DV-60 LED can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment(transmitters) and the models DV-60 LED are recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter/W.	Separation distance according to frequency of transmitter/m.		
	150 kHz~ 80 MHz $d = [\frac{3.5}{V_1}] \sqrt{p}$	80 MHz~800 MHz $d = [\frac{3.5}{E_1}] \sqrt{p}$	800 MHz~2.5 GHz $d = [\frac{7}{E_1}] \sqrt{p}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3

100	12	12	23
<p>For transmitter rated at a maximum output power not listed above, the recommended separation distance d in meters(m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts(W) given by transmitter' s manufacturer.</p> <p>Note1: At 80 MHz and 800 MHz, adopt the higher frequency range formula.</p> <p>Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

**The DV-60 LED light curing device has passed tests according IEC 60601-1-2, but it is no guarantee of immunity from electromagnetic interference. Avoid using The DV-60 LED light curing device in high electromagnetic environment.**

Manufacturer Name:URIT Medical Electronic Co., Ltd.

Address:No.D-07 Information Industry District,  
High-Tech Zone, Guilin, Guangxi 541004, P.R.China

Tel:+86(773)2288586 Fax:+86(773)2288560

Web:www.urit.com Email:service@uritest.com



Rev. 07/2023

Period of use:8 years

File No.:DVP050000096

Date of compilation:2023/07/31